




Aerolizer®

The Aerolizer works by placing a capsule into the device.

Loading a Dose of Medicine:

1. Take a capsule out of the wrapper.
2. Take the lid off the Aerolizer by twisting it.
3. Place the capsule in the hole in the center of the device.
Replace the lid.
4. Hold the Aerolizer with the mouthpiece up. Press the buttons on the sides of the Aerolizer once firmly. This action will pierce the capsule so you can inhale the medicine when you breathe in.



Breathing In the Medicine:

5. Before breathing in the medicine, turn your head away from the Aerolizer and breathe out.
6. Seal your lips around the Aerolizer mouthpiece.
7. Breathe in through your mouth. Remove the Aerolizer from your mouth.
8. Resume normal breathing.
9. Twist the lid of the Aerolizer to open. Throw away the used capsule. Always place the lid back on the Aerolizer when not in use to keep it clean and dry.